



TECHNIQUE CHART

MODEL NUMBER _____

SERIAL NUMBER _____

| X-RAY VIEW | mA | kVp | TIME | COMMENTS |
|------------|----|-----|------|----------|
| | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |

All X-Cel X-Ray stationary floor or table models have a fixed (28", 29" or 40") SID.

| <u>PULSES</u> | <u>SECONDS</u> | <u>PULSES</u> | <u>SECONDS</u> |
|---------------|----------------|---------------|----------------|
| 6 | .1 | 39 | .65 |
| 9 | .15 | 42 | .70 |
| 12 | .20 | 45 | .75 |
| 15 | .25 | 48 | .80 |
| 18 | .30 | 51 | .85 |
| 21 | .35 | 54 | .90 |
| 24 | .40 | 57 | .95 |
| 27 | .45 | 60 | 1.0 |
| 30 | .50 | 75 | 1.25 |
| 33 | .55 | 90 | 1.50 |
| 36 | .60 | 120 | 2.00 |